

# EMG-Controlled Exosuit Setup Guide

## 1 Project Overview

This guide provides complete setup instructions for the EMG-controlled exosuit system, which uses Myo armband data to control an exosuit via machine learning models.

## 2 Prerequisites

- Python 3.7+
- Android Studio
- Myo Armband
- Raspberry Pi (for deployment) or WSL (for development)

## 3 Step 1: Repository Setup

1. Clone the main repository:

```
1 git clone <https://gitlab.aibe.fau.de/airob/theses/research-lab-dawood-mughal>
2 cd <research-lab-dawood-mughal>
3
```

2. Install required Python dependencies:

```
1 pip install -r requirements.txt
2
```

## 4 Step 2: Training Server Setup

### 4.1 Training Server Configuration

1. Navigate to the TrainerUDPServer folder:

```
1 cd TrainerUDPServer
2
```

2. The file `wrist_exo_model_trainer.py` serves as the training server with UDP communication instead of traditional HTTP calls.
3. Key features of the training server:
  - Receives EMG recorded CSV files via UDP from Android
  - Trains machine learning models based on configurable feature sets
  - Implements `find_best_model()` function for model selection
  - Allows customization of feature sets within this function

## 4.2 Running the Training Server

1. Run the training server directly (not in WSL):

```
1 python wrist_exo_model_trainer.py
2
```

2. Configure firewall rules to allow ports: 3350, 3352, 3358, 12346, 12347
3. Note the IP address where the training server is running
4. Update the Android module with this IP address in `UdpMotorController.getTrainingServerIp()`




Figure 1: Training Server Configuration

## 5 Step 3: Motor Controller Setup

### 5.1 CANDLE and Motor Connections

1. Connect motors to each other in the designated sequence
2. Connect motors to the CANDLE device
3. Plug in the CANDLE device until the motor displays a blinking light indicating connection

### 5.2 Option 1: Raspberry Pi Deployment

1. Navigate to the controller directory:

```
1 cd pythonModule/controller
2
```

2. Modify configuration files:

- Set `TEST_MODE = False` in `exo_controller*.py`
- Set `TEST_MODE = False` in `main.py`
- Set `TEST_MODE = False` in `utils.py`

3. Run the main controller:

```
1 python main.py
2
```

## 5.3 Option 2: WSL Development Setup

### 1. Open WSL in VS Code:

- Open VS Code
- Press Ctrl + Shift + P and type "Remote-WSL: New Window"
- Select "Remote-WSL: New Window using WSL"
- In the new WSL window, open the path: /mnt/c/path/to/PythonModule/controller

### 2. Run main.py within WSL

### 3. Additionally run the UDP relay script outside WSL by making a new window as shown in Figure 2:

```
1 python udp_relay_for_wsl.py
2
```

### 4. Update WSL\_IP in udp\_relay\_for\_wsl.py with your WSL IP address

### 5. Note the system IP address (not WSL IP) for Android configuration




Figure 2: Motor Controller Setup

## 6 Step 4: Android Application Setup

### 6.1 Application Configuration

1. Open the project in Android Studio
2. Navigate to UdpMotorController in the utility folder
3. Update IP addresses in two locations:
  - getRaspiServerIp() - Motor controller IP
  - getTrainingServerIp() - Training server IP
4. Build and install the application on Android device

## 7 Step 5: System Operation

### 7.1 Initial Setup and Training

1. Launch the Android application
2. Accept required permissions
3. Scan for Myo Armband
4. Connect Myo armband and wear it properly
5. Click "Start Training Session" button




Figure 3: First Screen



Figure 4: Connect Myo Armband

## 7.2 Data Recording Process

1. Follow guided recording screen for EMG data collection
2. Perform movements as shown (flexion, extension, etc.)
3. After recording, choose Save data and train model:
  - Save data and train model
  - Save data only
  - Discard data
4. Select training algorithm (Ridge Regression or MLP)
5. Monitor training progress on screen




Figure 5: EMG Data Recording step




Figure 6: Model Selection

### 7.3 Select Trained Model and Activate

1. Return to home screen after completion
2. Use "Change/Select Model" button to choose trained model
3. Enable "Active Model" switch




Figure 7: Model Selection Interface

### 7.4 Exosuit Activation and Control

1. Click "Motor Settings" on main screen
2. Adjust motor parameters as needed
3. Save and apply settings (sends initial configuration)
4. Monitor connection status:
  - Initially: Disconnected
  - After settings applied: Status becomes READY TO START

5. Press "Start System" button to begin real-time control
6. Use activate model switch to pause sending predictions.
7. Use disconnect button to terminate connection




Figure 8: Motor Settings




Figure 9: Ready To Start System




Figure 10: Main Control Interface




Figure 11: Real-Time View

## 8 Troubleshooting

### 8.1 Common Issues

- **Firewall blocking ports:** Ensure ports 3350, 3352, 3358, 12346, 12347 are open. If there is a problem with UDP communication even after adding inbound rule try disabling Firewall for all.
- **Connection timeout:** Verify IP addresses in Android configuration
- **Training failures:** Check feature set configuration in `find_best_model()`
- **Motor not responding:** Confirm `TEST_MODE` is set to `False` in controller files

### 8.2 Connection Status Indicators

- **Red:** Disconnected
- **Yellow:** Connecting
- **Green:** Connected and active



## 9 Conclusion

This setup guide provides complete instructions for deploying the EMG-controlled exosuit system. Follow each step sequentially and refer to the accompanying figures for visual guidance. The system enables real-time control of exosuit movements based on EMG signal predictions from trained machine learning models.